

Analog Water Bath

- ❑ Easy to use, Economic Model
- ❑ Precise Temperature Control up to 99°C
- ❑ Built-in Safety and Convenience
- ❑ Three Models with Capacities from 6 to 22ℓ
- ❑ Seamless Stainless Steel Bath (SUS304)



Feature

Economic and Durable

general purpose water baths are ideal for industrial, pharmaceutical, clinical and biomedical use.

Precision Analog Controller

maintain water temperature from ambient +5°C to 99°C with $\pm 3.0^\circ\text{C}$ uniformity

Seamless Stainless Steel

corrosion resistant (SUS304) bath

Over Temperature Protection

ensure user safety

Perforated Sample Tray

protects heater and sensor from unexpected damages

Optional Accessories



Ring cover (LWB-R6)



Ring cover (LWB-R8)




Gabled Lid



Test Tube Rack

MODEL		LWB-206A	LWB-211A	LWB-222A
Dimensions (W x D x H mm)	Inner	300 x 155 x 150	300 x 240 x 150	500 x 295 x 150
	Outer	327 x 178 x 280	327 x 265 x 280	537 x 327 x 280
Capacity		6ℓ	11ℓ	22ℓ
Heater		500W	700W	1400W
Temperature	Range	Ambient + 5°C to 99°C		
	Set Accuracy	$\pm 0.5^\circ\text{C}$		
	Uniformity	$\pm 3.0^\circ\text{C}$		
Material	Interior	Seamless Stainless Steel (SUS304)		
	Exterior	Powder Coated Steel		
Safety		Over Temp. Cut-Off		
Electric Supply		110V 60Hz or 220V 50/60Hz		

 Custom made order is available
Specifications can be changed without prior notice for quality improvement.