

Basic Water Still



WD-2004F Basic Water Still /w Prefilter



STAINLESS STEEL EXTERIOR MODEL
LWD-104S / LWD-108S / LWD-204SF / LWD-208SF

Feature

- ❑ Space Saving Compact and Practical Design
- ❑ Maximum 7.5ℓ / hour Production Capacity
- ❑ Economical Energy and Water Consumption

Unique Design & Safety Feature

with hydraulic temperature controller and a water float switch provides automatic operation

Easy to Clean Up


by removing the condenser from the body

All Stainless Steel with Double Walled Chamber

finished with powder coated two colors

MODEL	WD-1004	WD-1008
	WD-2004F	WD-2008F
Dimensions (W x D x Hmm)	380 x 290 x 580	380 x 365 x 640
Capacity	3.5ℓ / hour	7.5ℓ / hour
Heater	3kW	6kW
Cooling Water	Approx. 30 ~ 50ℓ	Approx. 60 ~ 80ℓ
Controller	Hydraulic Thermostat	
Switch	Main Power with Pilot Lamp and Water Float Switch	
Sealing	Silicone Gasket	
Material	Inner	Stainless Steel
	Outer	Powder Coated Steel
	Heater	Stainless Steel Pipe
Safety	Low Water & Over Heating Cut-Out	
Electric Supply	220 V 50/60 Hz	

Oder Number	Descriptions	
WD-1004	Basic Water Still	3.5 ℓ / hour
WD-2004F	Basic Water Still /w Prefilter	3.5 ℓ / hour
WD-1008	Basic Water Still	7.5 ℓ / hour
WD-2008F	Basic Water Still /w Prefilter	7.5 ℓ / hour

 Custom made order is available
Specifications can be changed without prior notice for quality improvement.